



Taha Kimia Tajhiz Co.



Ofite Co.

Bulk Hardness Tester Datasheet

Drilling Fluids Testing/ Shear Strength Equipment



DRILLING FLUIDS EQUIPMENT

For over 30 years OFI Testing Equipment (OFITE) has provided instruments and reagents for testing drilling fluids, well cements, completion fluids, and wastewater. In addition to these product lines we also offer a range of instruments for core analysis. From our manufacturing facility in Houston, TX we provide customers all over the world with quality products and exceptional service.

Our drilling fluids product line includes innovative designs such as the Model 900 Viscosimeter, which showcases our ability to develop new technology to meet customer and industry demands. We also offer Retorts, Aging Cells, Roller Ovens, Mud Balances, Filter Presses, and all other instruments required to evaluate drilling fluid properties according to API Recommended Practice 13B-1 and 13B-2.

As an independent manufacturer and supplier, OFITE has one priority, our customers.

OFI TESTING EQUIPMENT, INC.
11302 Steeplecrest Dr.
Houston, TX 77065
877.837.8683
www.ofite.com
©Copyright OFITE 2016



Bulk Hardness Tester

The Bulk Hardness Tester (#150-87) is designed to evaluate the hardness of shale after exposure to fluids. The hardness of the shale can be related to the inhibitive properties of the fluid being evaluated. Shale that interacted with the fluids will become softer due to the adsorption of water, swelling, and dispersion of fine particles. This rock-fluid interaction can be linked to wellbore stability problems, including reduction in the compressive strength, spalling, or fracturing. In terms of integrity of drill cuttings, excessive softening and stickiness of the pieces of shale can produce mud rings in the annulus, sticking problems in the drilling assembly, and bit balling, among other problems.



Features

- Evaluates the hardness of shale samples after exposure to fluids
- Stainless steel construction
- Includes threaded piston, torque wrench, and Allen key